


UV-Vis Study of different fissure sealants



1st International Conference on Fluorescence based diagnostic of Oral Diseases—ICFOD, Montpellier 2012

Results

Material	ES	SML	UXT	DDS	DS
Filler content [%]	55	45	57	0.1	31.7
Transmission [%] 405 nm (emission)					
0.5 mm	54.2	20.6	14.5	73.8	27
1.0 mm	33.5	18	3.8	63	15.8
1.5 mm	19	16.1	0.8	53	8.2
2.0 mm	21	13.4	0.4	44.8	2
Transmission [%] 540 nm (irradiation)					
0.5 mm	82.2	33	26.1	85.6	51.4
1.0 mm	75	32.1	9.7	84.7	32.8
1.5 mm	61	23.8	3.1	81.6	17.1
2.0 mm	59.5	14.5	1.5	66.9	8.4
Transmission [%] 655 nm (emission, irradiation)					
0.5 mm	87	43.7	35.2	87.1	61.6
1.0 mm	83.3	37.1	15.7	86.7	44.7
1.5 mm	71.7	32.8	6	84.2	26.5
2.0 mm	72	18.4	3	70.4	17.1
Transmission [%] 675 nm (irradiation)					
0.5 mm	87.5	45.6	37	87.3	63
1.0 mm	84.2	37.6	17.1	87	46.5
1.5 mm	73	34.2	6.7	84.6	28.1
2.0 mm	73	19	3.4	70.8	18.8

Table 1: Filler content and UV-Vis data

PROPERTIES OF A HIGHLY FILLED CLEAR SEALANT FOR OPTICAL CARIES DIAGNOSIS

A. Barg¹, M. Danebrock¹, R. Maletz¹, A. Willner¹, F. Krause²
 1 VOCO GmbH, Cuxhaven, Germany
 2 Department of Operative and Preventive Dentistry, University of Bonn, Germany



Compressive Strength [MPa]




Figure 1: Compressive Strength

Flexural Strength [MPa]




Figure 2: Flexural Strength

Three Body Wear [µm]




Figure 3: Three body wear

Concerning
 als UXT, HI
 filled FV ar
 ES achieves
 the other m

tly supe-
 abrasion

nal high
 liagnosis
 iagnosis,
 a fissure

covered with a
 measure